

EC 40

DENTAL CEMENT REMOVER AQUEOUS ACID CLEANIG CONCENTRATE

Description

The concentrate removes fast and efficiently dental cements from bridges, crowns, dental and medical instruments made from: stainless-steel, ceramics, plastics, glass and rubber.

Removes dental cements including phosphate cements, metal oxides, lime and fluxing agents.

Not suitable for light metals (e.g. Al, Mg, Zn and corresponding alloys).

Application and dosage

Dosage: 5 – 10 % dilution (equivalent to 50 – 100 ml *EC 40* in 1 l solution with tap water). Clean in an ultrasonic bath for 5 – 10 min; then rinse thoroughly under running water.

Safety recommendations

EC 40 is classified as hazardous according to the regulation (EC) No 1272/2008 [GHS] (skin corrosive and serious eye damage, corrosive to metals [aluminium and light metals]). Observe also with respect to this the hints indicated in the Safety Data Sheet and always handle chemicals with care.

Physical-chemical characterisation

Density: 1182 g/l, pH (concentrate): 0.6 in original state.

Ingredients according to Annex VII, A, EC-Regulation 648/2004 (detergents): 15-30% phosphoric acid.

Disposal

The cleaning solution can be fed into the public sewage system after neutralisation; observe the local pH limit values and make sure that the contamination contained complies with the local sewage regulations.

Recommendation: use lime or soda lye for neutralisation.

European waste code: 20 01 14*, "acids".

Volumes, storage and transport

Available volumes: 1 litre PE-bottle; 10 litre and 25 litre HDPE-can

Do not store with alkalis. Provide acid-resistant floor. Store in closed original container at a temperature above -5°C, protected from heat and direct solar radiation. Store in a dry place.

Shelf life: 5 years from date of production (see stamp on label).

Classification for all means of transport: class 8, UN 1805.



Elma Schmidbauer GmbH

Gottlieb-Daimler-Str. 17, D-78224 Singen

e-mail: info@elma-ultrasonic.com

13.07.2016 / Version 2.0 / E:\PRODUKTINFO\englisch\word_GB\p_EC 40_GB.docx

Tel.: +49 7731 882-00

Fax: +49 7731 882-266

www.elma-ultrasonic.com