




Rotation inside your ultrasonic cleaning device.

Sieve cleaning module SRH 4/200



Rotating sieve cleaning

Angled positioned. Perfectly cleaned.



The SRH 4/200 offers an efficient yet gentle solution for cleaning analysis sieves. Combined with ultrasonic technology, it impresses with its sophisticated technology and visible results.

Benefits

- : Uniform cleaning of up to four analysis sieves (Ø up to 200/203 mm, 8") using ultrasonic technology and continuous rotation.
 - : Reliable cleaning and rinsing thanks to the angled sieve position.
 - : Cost-effective upgrade that can be retrofitted to an existing ultrasonic cleaning device without the need for tools (depending on the tank size).
-

Product information.

Order numbers

Sieve cleaning module	SRH 4 / 200
Order number 230V	1042706
Order number 115V	1044904


Mechanical data

Sieve cleaning module	Units	SRH 4 / 200
Max. ext. dimensions W / D / H	mm	675 / 295 / 330
Weight	kg	6.8

Modules per unit (maximum)

Sieve cleaning module	SRH 4 / 200
Elmasonic P / Select / Easy 300	1 x
Elmsonic Select 500	1 x
Elmsonic Select 900	2 x

Definition of dimensions





Elma Schmidbauer GmbH
Gottlieb-Daimler-Straße 17
78224 Singen (Germany)

www.elma-ultrasonic.com

