Product Information/Directions for use



elma tec clean A5

Description

elma tec clean A5 is an alkaline ultrasonic cleaner (powder) suitable for metal, artificial and natural rubber surfaces, gray cast iron and zinc die casting as well as for light metal alloys and plastic.

Removes scaled or resinous oil and grease, wax, residual lacquer and paint.

Check magnesium alloys and coatings for compatibility before use. Not suitable for plastics that are not alkali-proof or that are susceptible for stress cracking.

Application and Dosage

- Ultrasonic bath: Dosage: 3 5 % by weight Temperature: 50-80°C Cleaning time: 3-10 min.
- Warm up water in the tank to at least 50 °C, start the ultrasonic, add powder slowly while constantly stirring.
- After cleaning rinse the surfaces thoroughly.

Safety Recommendations

elma tec clean A5 is classified as hazardous according to the Regulation (EC) No 1272/2008 [GHS] (acutely toxic (inhalation) - category 4, skin corrosive and serious eye damage, respiratory irritant effect).

Observe also with respect to this the hints indicated in the safety data sheet and always handle chemicals with care.

Physical-chemical characterisation

- Bulk density: 0.92 g/ml pH (1 % solution): ~12.
- Ingredients according to Annex VII, A, Regulation (EC) No 648/2004 on detergents: 5-15 % anionic surfactants, <5 % non-ionic surfactants, >30 % phosphates.

Disposal

The surfactants in our product meet the criteria for biodegradation as laid down in Annex III of the Regulation (EC) No 648/2004 on detergents.

The used cleaning solution can be fed into the public sewage system after neutralisation; observe the local pH limit values and make sure that the contamination contained in the used cleaning solution complies with the local sewage regulations. For neutralisation we recommend: for 0.5 kg **elma tec clean A5** add step by step approx. 230 g water-free citric acid to the ultrasonic bath – do not use hydrochloric or sulphuric acid there!

European waste code: 20 01 29*, "detergents containing dangerous substances".

Volumes, Storage and Transport

- Available volume: 0.85 kg pot EFF 8000141 or 25 kg UN-approved cardboard drum EFF 8000140.
- Store in closed original container at a temperature between -5 °C and +30 °C, keep container dry and tightly closed.Do not store with acids.
- Shelf life: 5 years from date of production (see lot on label).
- Classification for all means of transport: class 8, UN 3253.

Accessories

■ Dosing cup: 250 ml vol./5 ml-scale ■ 8000643 • 1000 ml vol./10 ml-scale ■ 8000647.