

Elmasonic Accessory - Inserts


Holding clamps for Erlenmeyer flasks

- made of stainless steel and for mounting on the bottom of a cleaning basket
- the Erlenmeyer flask has a firm hold in the holding clamp due to the curved spring
- additional information on the number of Erlenmeyer flasks/shaking bottles with and without holding clamp in the Elma cleaning baskets → see FAQ "Erlenmeyer flasks"



	50 ml	250 ml	500 ml	1000 ml
Order number	1007553	1007555	1007557	1007558
List price	22 €	28 €	33 €	40 €
Ext. dimensions W/D/H (mm)	55/55/50	90/90/70	110/110/85	130/130/85

Dimension definition



All prices ex works, excluding packaging, plus VAT.

The drawings are merely included as a visualization aid and are not true to scale.

Technical specifications serve as orientation values and are subject to tolerances due to production and manufacturing. Proper operation can only be guaranteed when using original Elma accessories.