



Elmasteam Accessories - Medical Adapter Set

Set of different nozzle attachments

- connectable to an Elmasteam 8 med steam cleaner with standardized Luer-Lock connections
- allows the cleaning of components and instruments with difficult-to-access areas
- storage rack for 8 different attachments
- pliers for changing the attachments
- adapter Luer-Lock
- reliable cleaning for a wide range of applications

	Medical Adapter Set
Order number	1028820
List price	600,00 €
Ext. Dimensions W / D / H (mm)	130 / 70 / 560
Height rack H ₁ (mm)	40












All prices ex works, excluding packaging.

The drawings are meant as a visualization aid and are not true to scale.

Technical specifications serve as orientation values and are subject to tolerances due to production and manufacturing. Proper operation can only be guaranteed when using original Elma devices.

Product Information and Prices



	MIS hollow needle	MIS hollow needle with brush	Hollow needle attachment straight	Hollow needle attachment circumflex	Catheter attachment large	Catheter attachment medium	Hose connection large	Adapter Luer-Lock fem. / fem.	Pliers
Order number	1028822	1055030	1028823	1028824	1028825	1028826	1028827	1028828	1028831
List price	199,00 €	178,00 €	36,50 €	36,50 €	40,00 €	40,00 €	31,50 €	25,50 €	48,50 €
Total length L ₁ (mm)	525	338	103	88	44	44	26	21	115
Usable length L ₂ (mm)	501	315	84	61	29	29	12	9	16
Diameter D (mm)	2,5	2,5	2	2	2-9	2,5-6,5	3,5	9	--
Width W (mm)	--	--	--	--	--	--	--	--	29
Height H (mm)	--	--	--	--	--	--	--	--	15
									

All prices ex works, excluding packaging.

The drawings are meant as a visualization aid and are not true to scale.


Technical specifications serve as orientation values and are subject to tolerances due to production and manufacturing. Proper operation can only be guaranteed when using original Elma devices.

Elmasteam accessories - Medical adapter set

Dimension definition


MIS hollow needle

for tube shafts,
lateral steam outlet




MIS hollow needle

with brush, for tube shafts,
steam outlet behind brush hem




Hollow needle attachment

straight, front steam outlet



Hollow needle attachment

curved, front steam outlet




The drawings are meant as a visualization aid and are not true to scale.

Technical specifications serve as orientation values and are subject to tolerances due to production and manufacturing. Proper operation can only be guaranteed when using original Elma devices.


Elmasteam accessories - Medical adapter set

Dimension definition


Catheter attachment
large,
front steam outlet




Catheter attachment
medium,
front steam outlet




Hose connection
large,
front steam outlet



Adapter Luer-Lock
fem. / fem.,
front steam outlet



Pliers




The drawings are meant as a visualization aid and are not true to scale.

Technical specifications serve as orientation values and are subject to tolerances due to production and manufacturing. Proper operation can only be guaranteed when using original Elma devices.

Elmasteam accessories - Medical adapter set

Adapter for Elmasteam basic

Dimension definition



	Adapter
Order number	1091332
List price	168,00 €
Total length L_1 (mm)	33
Usable length L_2 (mm)	30
Diameter D (mm)	10
	

All prices ex works, excluding packaging.

The drawings are meant as a visualization aid and are not true to scale.

Technical specifications serve as orientation values and are subject to tolerances due to production and manufacturing. Proper operation can only be guaranteed when using original Elma devices.