Product Information and Prices





Elmasteam Accessories - Water canister

- Integrated level control
- Automatic filling of ElmaSteam 8 steam jet units with integrated pump
- Uninterrupted filling steam jet unit can remain pressurized

Order number	1070254
List price	237,00 €
Ext. Dimensions W / D / H (mm)	235 / 195 / 305
Capacity (I)	10
Hose length (m)	2

Dimension definition





All prices ex works, excluding packaging.

The drawings are meant as a visualization aid and are not true to scale.

Technical specifications serve as orientation values and are subject to tolerances due to production and manufacturing. Proper operation can only be guaranteed when using original Elma accessories.