

Product Information and Prices



Elmasteam Accessories - Elma Steam Descaler

- Descaler, citric acid-based, for Elmasteam 8
- 1 packet is easy and safe to dose
- Odorless and tasteless
- Removes and prevents the build-up of limescale and boiler scale
- Descaling extends the service life, saves energy and expensive repairs



Order numbers and list prices

Order number	8000878
List price	12,00 €

The prices indicated are net prices ex works, excluding transport costs, insurance, customs duties and taxes.

Product contents

- 5 packets of Elma Steam Descaler (45g each)
- Brief instructions descaling

Technical specifications serve as orientation values and are subject to tolerances due to production and manufacturing. Proper operation can only be guaranteed when using original Elma accessories.