


Elmasonic Table top units Accessories

Product Information and Prices



Contents

| | |
|---|----|
| - Modular basket system | 3 |
| Modular basket | 3 |
| Level set | 5 |
| Dividing wall set | 7 |
| Holder set | 9 |
| Silicone holder set | 11 |
| Elma starter Set | 12 |
| Elma Pin Set | 13 |
| Silicone holder set (Bottom comb / Side comb) | 15 |
| Basket for small parts | 17 |
| - Baskets | 19 |
| Basket | 19 |
| Stacking basket | 21 |
| - Covers | 22 |
| Stainless steel lid | 22 |
| Plastic lid | 24 |
| Noise protection covers | 26 |
| - Holders and product carriers | 27 |
| Equipment holder | 27 |
| Test tube holder | 28 |
| Holder for reaction tubes | 30 |
| Hooks holder | 32 |
| Ring stand | 34 |
| Holding clamps for Erlenmeyer flasks | 35 |
| Dental-Tray holder | 38 |
| - Inserts | 40 |
| Acid resistant insert tray | 40 |
| BD-Steribox | 42 |
| Cleaning cup | 43 |
| Cleaning glasses | 45 |
| Immersion baskets | 53 |
| Hostal tub | 55 |
| Silicone mats | 57 |
| Small parts holder | 59 |
| - Practicals | 60 |
| Cooling coil | 60 |
| Drain sieve | 61 |
| Sieve cleaning module SRH 4 / 200 | 62 |
| - Elma Mutebox | 64 |
| - Elma Sink | 67 |
| - Elma Dry | 70 |


Product Information and Prices



Elmasonic Accessories - Modular basket system

Modular basket

- stainless steel with plastic handles (10 - 120 & 180 - 300) or Rilsan coating (150)
- ergonomic handles
- autoclavable
- Bayonet holes in the bottom for mounting holders (e.g. Elma Pin, silicone holder bottom comb)
- Square openings in the side walls for mounting silicone holders



| Basket size | 10 | 20 | | 30 | 40 | 60 | 80 | 100 | 120 | 150 | 180 | 300 | | | | | | |
|---|---------------|----------------|----------------|----------------|----------------|----------------|----------------|-----------------|----------------|-----------------|-----------------|-----------------|--|--|--|--|--|--|
| Order number | 1117923 | 1119922 | | 1108916 | 1113045 | 1113050 | 1111838 | 1113387 | 1113390 | 1113110 | 1113404 | 1113072 | | | | | | |
| List price | 45,00 € | 47,50 € | | 50,00 € | 58,00 € | 70,00 € | 130,00 € | 150,00 € | 150,00 € | 195,00 € | 175,00 € | 195,00 € | | | | | | |
| Int. dim. of basket W / D / H (mm) | 155 / 70 / 35 | 110 / 105 / 50 | 195 / 105 / 50 | 195 / 105 / 95 | 255 / 105 / 95 | 455 / 105 / 55 | 250 / 195 / 90 | 250 / 195 / 135 | 450 / 255 / 55 | 280 / 255 / 135 | 450 / 255 / 135 | 450 / 255 / 135 | | | | | | |
| Chamfer V x W (mm) | -- | 20 x 20 | | | | | 25 x 25 | | | | | | | | | | | |
| Hole diameter bottom Y x Z (mm) | 11 x 15 | | | | | | | | | | | | | | | | | |
| Hole diameter sides Y ₁ x Z ₁ (mm) | -- | 10 x 10 | | | | | | | | | | | | | | | | |
| Loading weight max. (kg) | 1,0 | 1,0 | | 1,0 | 3,0 | 5,0 | 5,0 | 6,0 | 7,0 | 7,0 | 8,0 | 10,0 | | | | | | |
| Suitable for unit / bath size | 10 | 15 | 20 | 30 | 40 | 60 | 80 | 100 | 120 | 150 | 180 | 300 | | | | | | |
| • Elmasonic Easy | • | • | • | • | • | • | • | • | • | • | • | • | | | | | | |
| • Elmasonic Med | | | | • | | • | | • | • | | • | • | | | | | | |
| • Elmasonic Select | | | | • | • | • | | • | • | • | • | • | | | | | | |
| • Elmasonic S | • | • | | • | • | • | | • | • | • | • | • | | | | | | |
| • Elmasonic P | | | | • | | • | • | | • | | • | • | | | | | | |
| • Elmady | | | | • | | | | | • | | | • | | | | | | |
| • Elma Sink | | | | | • | • | | | | | • | • | | | | | | |



All prices ex works, excluding packaging.

Technical specifications serve as orientation values and are subject to tolerances due to production and manufacturing. Proper operation can only be guaranteed when using original Elma devices.

Elmasonic Accessories - Modular basket system

Modular basket

Dimension definition



The drawings are meant as a visualization aid and are not true to scale.

Technical specifications serve as orientation values and are subject to tolerances due to production and manufacturing. Proper operation can only be guaranteed when using original Elma devices.

Elmasonic Accessories - Modular basket system

Level set

- Level made of stainless steel
- for horizontal division of the modular basket
- Product contents Level set:
 - Level size S and M: 1x Level + 5x Silicone holder per set
 - Level size L and XL: 1x Level + 10x Silicone holder per set
- maximum one horizontal division per modular basket allowed
- the level can be used at different heights
- autoclavable
- Bayonet holes in the bottom for mounting holders (e.g. Elma Pin, silicone holder bottom comb)

| Level size | S | M | L | XL |
|--|---------------|---------------|---------------|---------------|
| Order number | 1125225 | 1125226 | 1125227 | 1125228 |
| List price | 15,00 € | 15,00 € | 30,00 € | 35,00 € |
| Ext. dimension W / D / H (mm) | 195 / 105 / 1 | 255 / 105 / 1 | 250 / 195 / 1 | 285 / 250 / 1 |
| Hole diameter Dia (mm) | 11 x 15 | | | |
| Weight per level (g) | 110 | 150 | 280 | 410 |
| Max. loading weight per basket size and level (kg) | 1,0 | 1,0 | 2,0 | 2,0 |

All prices ex works, excluding packaging.

Technical specifications serve as orientation values and are subject to tolerances due to production and manufacturing. Proper operation can only be guaranteed when using original Elma devices.




| Level size | S | M | L | XL |
|---------------------------|---------------------|----------|----------|---------------------|
| Suitable for basket size: | | | | |
| 30 | 1 set | | | |
| 40 | 1 set | | | |
| 60 | | | 1 set | |
| 80 | 2 sets side by side | | | |
| 100 | | | | 1 set |
| 120 | | | | 1 set |
| 150 | | | | 2 sets side by side |
| 180 | | | | 1 set |
| 300 | | | | 2 sets side by side |

Elmasonic Accessories - Modular basket system


Level set

Dimension definition




Example of application

Basket size 60 with Level set size M



Basket size 300 with 2 Level sets size L side by side




The drawings are meant as a visualization aid and are not true to scale.

Technical specifications serve as orientation values and are subject to tolerances due to production and manufacturing. Proper operation can only be guaranteed when using original Elma devices.

Elmasonic Accessories - Modular basket system

Dividing wall set

- Dividing wall made of stainless steel
- for vertical division of the modular basket
- Product contents Dividing wall set: 2x Dividing wall + 10x Silicone holder per set
- autoclavable
- Square openings in the side walls for mounting silicone holders



| Dividing wall size | S low | S high | M | L low | L high |
|----------------------------------|--------------|---------------|--------------|--------------|---------------|
| Order number | 1125220 | 1125221 | 1125223 | 1125222 | 1125224 |
| List price | 25,00 € | 40,00 € | 35,00 € | 35,00 € | 40,00 € |
| Ext. dimension W / D / H (mm) | 105 / 50 / 1 | 105 / 90 / 1 | 200 / 60 / 1 | 200 / 90 / 1 | 255 / 130 / 1 |
| Hole diameter Dia (mm) | 10 x 10 | | | | |
| Weight per dividing wall (g) | 70 | 40 | 130 | 100 | 230 |

| Dividing wall size | S low | S high | M | L low | L high |
|---------------------------|--------------|---------------|----------|--------------|---------------|
| Suitable for basket size: | | | | | |
| 20 | 1 Set | | | | |
| 30 | 2 Sets | | | | |
| 40 | 2 Sets | 2 Sets | | | |
| 60 | 2 Sets | 2 Sets | | | |
| 80 | 4 Sets | 4 Sets | | | |
| 100 | 2 Sets | 2 Sets | 2 Sets | | |
| 120 | 2 Sets | 2 Sets | 2 Sets | | |
| 150 | 4 Sets | | | 4 Sets | |
| 180 | 2 Sets | 2 Sets | | | 2 Sets |
| 300 | 4 Sets | 4 Sets | | | 2 Sets |


All prices ex works, excluding packaging.

Technical specifications serve as orientation values and are subject to tolerances due to production and manufacturing. Proper operation can only be guaranteed when using original Elma devices.

Elmasonic Accessories - Modular basket system


Dividing wall set

Dimension definition




Example of application

Basket size 300 with 2x Dividing wall set L high + 1x Dividing wall set S high



Basket size 30 with 2x Dividing wall set S low




The drawings are meant as a visualization aid and are not true to scale.

Technical specifications serve as orientation values and are subject to tolerances due to production and manufacturing. Proper operation can only be guaranteed when using original Elma devices.

Elmasonic Accessories - Modular basket system

Holder set


- made of stainless steel
- Holder hooks Rilsan-coated
- Holder hooks with 10 hooks
- Holder hooks for holding e.g. jewelry and watch parts
- Holder shaft with 4 shaft
- Holder shaft for holding e.g. Impression tray
- Product contents Holder set: 1x Holder shaft, 1x Holder hooks, 5x Silicone holder per set




| | Holder set |
|---|---------------|
| Order number | 1125219 |
| List price | 40,00 € |
| Ext. dimension Holder shaft W / W / H (mm) | 105 / 3 / 15 |
| Ext. dimension Holder hooks W ₁ / D ₁ / H ₁ (mm) | 110 / 55 / 25 |
| Weight bracket shaft (g) | 7 |
| Weight holder hooks (g) | 9 |

dimension definition

Holder shaft



Holder hooks



All prices ex works, excluding packaging.

The drawings are meant as a visualization aid and are not true to scale.


Technical specifications serve as orientation values and are subject to tolerances due to production and manufacturing. Proper operation can only be guaranteed when using original Elma Units.

Elmasonic Accessories - Modular basket system

Holder set

Example of application

Basket size 60 with 1x Holder set



The drawings are meant as a visualization aid and are not true to scale.

Technical specifications serve as orientation values and are subject to tolerances due to production and manufacturing. Proper operation can only be guaranteed when using original Elma devices.

Elmasonic Accessories - Modular basket system



Silicone holder set

- for mounting levels, dividing walls and supports
- Product contents Silicone holder set: 10x Silicone holder per set
- easy handling
- easy positioning



dimension definition

| | Silicone holder set |
|-----------------------------------|---------------------|
| Order number | 1125208 |
| List price | 15,00 € |
| Ext. dimension Dia / H (mm) | 20 / 20 |
| Weight per silicone holder (g) | 4 |



All prices ex works, excluding packaging.

The drawings are meant as a visualization aid and are not true to scale.


Technical specifications serve as orientation values and are subject to tolerances due to production and manufacturing. Proper operation can only be guaranteed when using original Elma Units.

Elmasonic Accessories - Modular basket system

Elma starter Set

- Product contents Elma starter Set: 3x each Elma Pin 0741 and 0772 per set
2x each Elma Pin 0509, 0158, 0767, 0613+0614, 0021R and 0001R per set
2x silicone holder side comb and bottom comb each per set
- for further features on the individual Elma pins and silicone holders, see page 9 onwards

| | Elma starter Set |
|--------------|------------------|
| Order number | 1124042 |
| List price | 50,00 € |



All prices ex works, excluding packaging.

The drawings are meant as a visualization aid and are not true to scale.

Technical specifications serve as orientation values and are subject to tolerances due to production and manufacturing. Proper operation can only be guaranteed when using original Elma Units.

Product Information and Prices



Elmasonic Accessories - Modular basket system

Elma Pin Set

- made of plastic
- Color: Black
- For mounting and holding various components
- Product contents Elma Pin Set: 10x Elma Pin per set
- Bayonet lock for easy mounting on the bottom of the baskets and on the levels



| Elma Pin | 0509 | 0158 | 0767 | 0613 + 0614 | 0741 | 0772 | 0021R | 0001R |
|--|---------|---------|---------|--------------|---------|---------|---------|---------|
| Order number | 1125209 | 1125210 | 1125211 | 1125212 | 1125214 | 1125213 | 1125216 | 1125215 |
| List price | 15,00 € | 15,00 € | 15,00 € | 50,00 € | 15,00 € | 15,00 € | 15,00 € | 15,00 € |
| Ext. dimension Dia / H (mm) | 30 / 35 | 50 / 40 | 85 / 65 | 55 / 45 | 75 / 45 | 25 / 30 | 25 / 55 | 40 / 75 |
| int. dimension Dia ₁ / H ₁ (mm) | 20 / 20 | 40 / 25 | 50 / 35 | 0,35 / 25 | 25 / 25 | 4 / 15 | 5 / 40 | 15 / 40 |
| Outer overlay Dia ₂ (mm) | -- | -- | -- | 15 | 40 | -- | -- | -- |
| Weight per Elma Pin (g) | 3 | 6 | 7 | 6 | 2 | 3 | 4 | 8 |



All prices ex works, excluding packaging.



Technical specifications serve as orientation values and are subject to tolerances due to production and manufacturing. Proper operation can only be guaranteed when using original Elma devices.

Elmasonic Accessories - Modular basket system



Elma Pin Set

Dimension definition



Elma Pin 0509 & 0158





Elma Pin 0767





Elma Pin 0613 & 0614



Elma Pin 0741 & 0772



Elma Pin 0021R & 0001R



The drawings are meant as a visualization aid and are not true to scale.

Technical specifications serve as orientation values and are subject to tolerances due to production and manufacturing. Proper operation can only be guaranteed when using original Elma Units.



Elmasonic Accessories - Modular basket system

Silicone holder set (Bottom comb / Side comb)

- for holding various instruments, containers (e.g. flasks) and workpieces
- Scope of delivery silicone holder set bottom comb: 10x silicone holder bottom comb each per set
- Scope of delivery silicone holder set side comb: 5x silicone holder side comb each per set
- Color: White



dimension definition



| | Bottom comb | Side comb |
|--------------------------------|---------------|---------------|
| Order number | 1125217 | 1125218 |
| List price | 60,00 € | 60,00 € |
| Ext. dimension W / D / H (mm) | 105 / 10 / 35 | 150 / 10 / 45 |
| Weight per silicone holder (g) | 35 | 71 |

All prices ex works, excluding packaging.

The drawings are meant as a visualization aid and are not true to scale.


Technical specifications serve as orientation values and are subject to tolerances due to production and manufacturing. Proper operation can only be guaranteed when using original Elma Units.

Elmasonic Accessories - Modular basket system


Silicone holder set (Bottom comb / Side comb)

Example of application

Basket size 150 with 2x Bottom comb



Basket size 60 with 1x Side comb



The drawings are meant as a visualization aid and are not true to scale.

Technical specifications serve as orientation values and are subject to tolerances due to production and manufacturing. Proper operation can only be guaranteed when using original Elma Units.

Product Information and Prices



Elmasonic Accessories - Modular basket system

Basket for small parts

- Basket for cleaning small parts
- made of stainless steel
- compatible with Elma pins



| | Basket 1-piece 1 / 3 height | Basket 5-pieces 1 / 3 height | Miniature basket | Cover for Miniature basket | Cover for basket | Basket 4-pieces 1 / 3 height | Basket 9-pieces 1 / 3 height | Basket 1-piece 1 / 1 height |
|---|---------------------------------------|--|-------------------------|--------------------------------------|----------------------------|--|--|---------------------------------------|
| Order number | 3000006012 | 3000006014 | 3000006017 | 3000006016 | 3000006018 | 3000006013 | 3000006015 | 3000009104 |
| List price | 21,00 € | 40,00 € | 13,50 € | 17,00 € | 13,50 € | 35,00 € | 62,00 € | 40,00 € |
| Ext. dimensions Dia / H (mm) | 64 / 16 | 64 / 16 | 21 / 12 | 20 / 4 | 64 / 2 | 64 / 16 | 64 / 16 | 64 / 42 |
| Int. dimensions Dia ₁ / H ₁ (mm) | 60 / 15 | 60 / 15 | 18 / 11 | -- | -- | 60 / 15 | 60 / 15 | 60 / 41 |

All prices ex works, excluding packaging.

The drawings are meant as a visualization aid and are not true to scale.


Technical specifications serve as orientation values and are subject to tolerances due to production and manufacturing. Proper operation can only be guaranteed when using original Elma devices.

Elmasonic Accessories - Modular basket system


Basket for small parts

Dimension definition

Baskets



Covers



The drawings are meant as a visualization aid and are not true to scale.

Technical specifications serve as orientation values and are subject to tolerances due to production and manufacturing. Proper operation can only be guaranteed when using original Elma devices.


Product Information and Prices



Elmasonic Accessories - Baskets

Basket

- made of stainless steel and with Rilsan-coated handles
- the basket 200 (order no. 1108458) also serves as a holder for sterilisation baskets
 - DIN 1 / 1 short (W / D - 485 / 253 mm) (H - max. 100 mm)
 - DIN 3 / 4 (W / D - 406 / 253 mm) (H - max. 100 mm)



| | 70 | 100 medium | 120 medium | 100 flat | 120 flat | 200 | 500 | 900 |
|--|----------------|----------------------|----------------------|--------------------|--------------------|----------------|-----------------|-----------------|
| Order number | 1004247 | 1071801 | | 1077166 | | 1108458 | 1107675 | 1107641 |
| List price | 168,00 € | 190,00 € | | 165,00 € | | 195,00 € | 305,00 € | 405,00 € |
| Basket int. dim. W / D / H (mm) | 465 / 105 / 50 | 255 / 200 / 45 | | 255 / 200 / 25 | | 490 / 255 / 95 | 460 / 300 / 195 | 545 / 460 / 190 |
| Mesh size Y x Z (mm) / Hole diameter Dia (mm) | | 9 x 1 | | | | 11 | | 11 |
| Loading weight max. (kg) | 5,0 | | 6,0 | | | 4,0 | 20,0 | 30,0 |
| Suitable for | 70 | 100 | 120 | 100 | 120 | 200 | 500 | 900 |
| • Elmasonic Easy | | • | • | • | • | | | |
| • Elmasonic Med | | • | • | • | • | • | • | • |
| • Elmasonic Select | | • | • | • | • | | • | • |
| • Elmasonic S | | • | • | • | • | | • | • |
| • Elmasonic P | • | | • | | • | | | |
| • Elmady TD | | | • | | • | | | |

All prices ex works, excluding packaging.


Technical specifications serve as orientation values and are subject to tolerances due to production and manufacturing. Proper operation can only be guaranteed when using original Elma devices.

Elmasonic Accessories - Baskets

Basket


Dimension definition

Basket internal dimensions





Hole diameter

Dia = Hole diameter (basket sizes 200, 500, 900)



Mesh size

Y = Mesh size; Z = Wire diameter (basket sizes 70, 100 medium / flat,
120 medium / flat)



The drawings are meant as a visualization aid and are not true to scale.

Technical specifications serve as orientation values and are subject to tolerances due to production and manufacturing. Proper operation can only be guaranteed when using original Elma devices.

Elmasonic Accessories - Baskets

Stacking basket

- made of stainless steel and with Rilsan-coated handles



| | 180 |
|--|----------------|
| Order number | 2070880000 |
| List price | 405,00 € |
| Basket int. dim. lower basket W / D / H (mm) | 270 / 240 / 50 |
| Basket int. dim. upper basket W / D / H (mm) | 265 / 235 / 50 |
| Mesh size Y x Z (mm) | 9 x 1 |
| Loading weight max. (kg) | 8,0 |

All prices ex works, excluding packaging.

Technical specifications serve as orientation values and are subject to tolerances due to production and manufacturing. Proper operation can only be guaranteed when using original Elma devices.

Elmasonic Accessories - Covers

Stainless steel lid

- reduce the noise of the ultrasonic cleaning
- save energy and costs by reducing waste heat and collecting / recirculating the water vapour



| | 30 | 40 | 60 | 100 | 120 | 150 | 180 | 300 | 500 | 900 | | | | | |
|---------------------------------------|----------------|-----------|----------------|------------|----------------|------------|----------------|------------|----------------|------------|----------------|--|----------------|--|----------------|
| Order number | 1111930 | | 1111931 | | 1111932 | | 1111933 | | 1111934 | | 1111933 | | 1082096 | | 1082091 |
| List price | 50,00 € | | 65,00 € | | 86,00 € | | 160,00 € | | 108,00 € | | 160,00 € | | 101,00 € | | 155,00 € |
| External dimensions W / D / H (mm) | 265 / 160 / 55 | | 325 / 175 / 55 | | 325 / 265 / 55 | | 530 / 325 / 55 | | 355 / 325 / 55 | | 530 / 325 / 55 | | 540 / 370 / 50 | | 650 / 550 / 50 |
| Suitable for | 30 | 40 | 60 | 100 | 120 | 150 | 180 | 300 | 500 | 900 | | | | | |
| • Elmasonic Easy | • | • | • | • | • | | • | • | | | | | | | |
| • Elmasonic Med | • | | • | • | • | | • | • | • | • | | | | | |
| • Elmasonic Select | • | • | • | • | • | • | • | • | • | • | | | | | |
| • Elmasonic S | • | • | • | • | • | • | • | • | • | • | | | | | |
| • Elmasonic P | • | | • | • | • | | • | • | | | | | | | |



All prices ex works, excluding packaging.

Technical specifications serve as orientation values and are subject to tolerances due to production and manufacturing. Proper operation can only be guaranteed when using original Elma accessories.

Elmasonic Accessories - Covers

Stainless steel lid

Dimension definition



The drawings are meant as a visualization aid and are not true to scale.

Technical specifications serve as orientation values and are subject to tolerances due to production and manufacturing. Proper operation can only be guaranteed when using original Elma devices.

Product Information and Prices



Elmasonic Accessories - Covers

Plastic lid

- available in cobalt blue (RAL 5013)*¹, ultramarine blue (RAL 5002)*² and silver grey (RAL 7001)*³
- reduce the noise of the ultrasonic cleaning
- save energy and costs by collecting and recirculating evaporating water from the ultrasonic bath
- serve upside down as a drip tray for the basket and cleaned parts



| | 10 | 15 | 20 | 30 | 40 | 60 | 70 | 80 | 100 | 120 | 150 | 300 | 180 | 200 | | |
|--|----------------|-----------|----------------|-----------|----------------|-----------|----------------|-----------|----------------|------------|----------------|------------|----------------|------------|----------------|----------------|
| Order number cobalt blue (RAL 5013)* ¹ | | | | 1101126 | | 1101129 | | 1110200 | | 1101133 | | 1101135 | | 1101141 | | 1106428 |
| Order number ultramarine blue (RAL 5002)* ² | 1003278 | | 1003280 | | 1003281 | | 1003282 | | 1003283 | | 1003285 | | 1003287 | | 1003286 | -- |
| Order number silver grey (RAL 7001)* ³ | 1071693 | | 1071694 | | 1071695 | | 1071696 | | 1110202 | | 1071697 | | 1071699 | | 1071698 | -- |
| List price | 14,50 € | | 20,00 € | | 23,00 € | | 28,00 € | | 37,00 € | | 37,00 € | | 45,50 € | | 39,00 € | 52,00 € |
| Ext. dimensions W / D / H (mm) | 235 / 115 / 20 | | 220 / 175 / 20 | | 305 / 175 / 20 | | 270 / 190 / 20 | | 570 / 170 / 20 | | 380 / 280 / 20 | | 575 / 335 / 20 | | 405 / 340 / 20 | 600 / 365 / 20 |
| Suitable for | 10 | 15 | 20 | 30 | 40 | 60 | 70 | 80 | 100 | 120 | 150 | 300 | 180 | 200 | | |
| • Elmasonic Easy | • | | • | • | | • | | | • | • | | • | | • | | |
| • Elmasonic Med | | | | • | | • | | | • | • | | • | | • | • | |
| • Elmasonic Select | | | | • | • | • | | • | • | • | • | • | | • | | |
| • Elmasonic S | • | • | | • | • | • | | • | • | • | • | • | | • | | |
| • Elmasonic P | | | | • | | • | • | | • | • | | • | | • | | |


All prices ex works, excluding packaging.

Technical specifications serve as orientation values and are subject to tolerances due to production and manufacturing. Proper operation can only be guaranteed when using original Elma devices.

Elmasonic Accessories - Covers

Plastic lid

Dimension definition



The drawings are merely included as a visualization aid and are not true to scale.

Technical specifications serve as orientation values and are subject to tolerances due to production and manufacturing. Proper operation can only be guaranteed when using original Elma devices.

Elmasonic Accessory - Covers



Noise protection covers

- reduce the noise of the ultrasonic cleaning by 40 % (500) and 54% (900)* and retain process vapours
- save energy and costs by reducing waste heat and collecting / recirculating the water vapour
- available as a retrofit kit for self-assembly and space-savingly attached to the ultrasonic device



| | 500 | 900 |
|---|-----------------|-----------------|
| Order number | 1100349 | 1100292 |
| List price | 455,00 € | 670,00 € |
| External dimensions W / D / H (mm) | 580 / 495 / 180 | 680 / 665 / 180 |
| Height installed H ₁ (mm) | 65 | |
| Suitable for | 500 | 900 |
| • Elmasonic Med | • | • |
| • Elmasonic Select | • | • |

Dimension definition



*Only when closed and without basket / load.

All prices ex works, excluding packaging.

The drawings are merely included as a visualization aid and are not true to scale.

Technical specifications serve as orientation values and are subject to tolerances due to production and manufacturing. Proper operation can only be guaranteed when using original Elma accessories.

Elmasonic Accessories - Holders and product carriers


Equipment holder

- made of stainless steel and with Rilsan-coated handles
- for holding laboratory vessels e.g. glass flasks up to 1000ml



| Order number | 1049704 |
|---|---------------|
| List price | 66,00 € |
| Basket int. dimension W / D / H (mm) | 195 / 95 / 75 |
| Suitable for | |
| Elmasonic Med, Select, Easy, P | |
| 60 | 1 x |
| 120 | 2 x |
| 180 | 2 x |
| 300 | 4 x |

Dimension definition



All prices ex works, excluding packaging.

The drawings are meant as a visualization aid and are not true to scale.

Technical specifications serve as orientation values and are subject to tolerances due to production and manufacturing. Proper operation can only be guaranteed when using original Elma Units.

Product Information and Prices



Elmasonic Accessories - Holders and product carriers

Test tube holder

- made of stainless steel
- for test tubes Ø 16 mm / Ø 25 mm



| | 30 | 40 | 60 | 80 | 100 | 120 | 150 | 180 | 200 | 300 |
|---------------------------------------|-----------|---------------|-----------|-----------|---------------|------------|---------------------------|---------------|------------|------------|
| Order number | | 1111973 | | | 1112084 | | | 1112123 | | |
| List price | | 65,00 € | | | 86,00 € | | | 100,00 € | | |
| External dimensions W / D / H (mm) | | 200 / 90 / 95 | | | 300 / 90 / 95 | | | 360 / 90 / 95 | | |
| Height hangers H ₁ (mm) | | 70 | | | | | 70 | | | |
| Diameter hole Dia (mm) | | 17 | | | | | 17 / 26 | | | |
| Max. number of test tubes | | 16 x Ø 16 mm | | | | | 8 x Ø 16 mm + 7 x Ø 25 mm | | | |
| Suitable for | 30 | 40 | 60 | 80 | 100 | 120 | 150 | 180 | 200 | 300 |
| • Elmasonic Easy | 1 x | 1 x | 2 x | 4 x | 2 x | 2 x | | 2 x | | 4 x |
| • Elmasonic Med | 1 x | | 2 x | 4 x | 2 x | 2 x | | 2 x | 4 x | 4 x |
| • Elmasonic Select | 1 x | 1 x | 2 x | 4 x | 2 x | 2 x | 4 x | 2 x | | 4 x |
| • Elmasonic S | 1 x | 1 x | 2 x | 4 x | 2 x | 2 x | 4 x | 2 x | | 4 x |
| • Elmasonic P | 1 x | | 2 x | 4 x | | 2 x | | 2 x | | 4 x |


All prices ex works, excluding packaging.

Technical specifications serve as orientation values and are subject to tolerances due to production and manufacturing. Proper operation can only be guaranteed when using original Elma accessories.

Elmasonic Accessories - Holders and product carriers

Test tube holder

Dimension definition



Example of application



The drawings are meant as a visualization aid and are not true to scale.

The parts shown are for illustration purposes only and are not included in the scope of delivery.

Technical specifications serve as orientation values and are subject to tolerances due to production and manufacturing. Proper operation can only be guaranteed when using original Elma accessories.

Elmasonic Accessories - Holders and product carriers

Holder for reaction tubes

- made of stainless steel
- for reaction tubes 1,5 oder 2 ml



| | 100 | 120 | 150 | 180 | 200 | 300 |
|---------------------------------------|---------------|------------|------------|---------------|------------|------------|
| Order number | 1112158 | | | 1112152 | | |
| List price | 80,00 € | | | 100,00 € | | |
| External dimensions W / D / H (mm) | 300 / 90 / 70 | | | 360 / 90 / 70 | | |
| Height hangers H ₁ (mm) | | | 35 | | | |
| Diameter hole Dia (mm) | | | 11 | | | |
| Suitable for | 100 | 120 | 150 | 180 | 200 | 300 |
| • Elmasonic Easy | 2 x | 2 x | | 2 x | | 4 x |
| • Elmasonic Med | 2 x | 2 x | | 2 x | 4 x | 4 x |
| • Elmasonic Select | 2 x | 2 x | 4 x | 2 x | | 4 x |
| • Elmasonic S | 2 x | 2 x | 4 x | 2 x | | 4 x |
| • Elmasonic P | | 2 x | | 2 x | | 4 x |


All prices ex works, excluding packaging.

Technical specifications serve as orientation values and are subject to tolerances due to production and manufacturing. Proper operation can only be guaranteed when using original Elma accessories.

Elmasonic Accessories - Holders and product carriers

Holder for reaction tubes

Dimension definition



Example of application



The drawings are meant as a visualization aid and are not true to scale.

The parts shown are for illustration purposes only and are not included in the scope of delivery.

Technical specifications serve as orientation values and are subject to tolerances due to production and manufacturing. Proper operation can only be guaranteed when using original Elma accessories.

Elmasonic Accessories - Holders and Product Carriers

Hooks holder

- with 10 hooks that can be flexibly attached
- for cleaning e.g. rings, chains, etc.
- plastic support surfaces to prevent scratching of the tray and the goods



| Unit size | 60 | 100 | 120 | 150 | 180 | 300 |
|---|-----|-----|-----|---------------|-----|-----|
| Order number | | | | 1117432 | | |
| List price | | | | 22,00 € | | |
| Ext. dimension W / D / H (mm) | | | | 355 / 60 / 75 | | |
| Support height H ₁ (mm) | | | | 55 | | |
| Dimensions hooks Dia / D ₁ (mm) | | | | 1,7 / 55 | | |
| Max. length with 10 hooks W ₁ (mm) | | | | 250 | | |
| Suitable for | 60 | 100 | 120 | 150 | 180 | 300 |
| • Elmasonic Easy | 1 x | 2 x | 2 x | 5 x | 4 x | 5 x |
| • Elmasonic Med | 1 x | 2 x | 2 x | | 4 x | 5 x |
| • Elmasonic Select | 1 x | 2 x | 2 x | 5 x | 4 x | 5 x |
| • Elmasonic S | 1 x | 2 x | 2 x | 5 x | 4 x | 5 x |
| • Elmasonic P | 1 x | | 2 x | | 4 x | 5 x |



All prices ex works, excluding packaging.

Technical specifications serve as orientation values and are subject to tolerances due to production and manufacturing. Proper operation can only be guaranteed when using original Elma accessories.

Elmasonic Accessories - Holders and Product Carriers

Hooks holder

Dimension definition



The drawings are meant as a visualization aid and are not true to scale.

The parts shown are for illustration purposes only and are not included in the scope of delivery.

Technical specifications serve as orientation values and are subject to tolerances due to production and manufacturing. Proper operation can only be guaranteed when using original Elma devices.

Elmasonic Accessories - Holders and Product Carriers

Ring stand

- made of stainless steel
- with 8 holding hooks



| | |
|-------------------------|------------|
| | |
| Order number | 2070460000 |
| List price | 29,00 € |
| Dimensions Dia / H (mm) | 80 / 160 |

Dimension definition



All prices ex works, excluding packaging.

The drawings are meant as a visualization aid and are not true to scale.

Technical specifications serve as orientation values and are subject to tolerances due to production and manufacturing. Proper operation can only be guaranteed when using original Elma devices.

Elmasonic Accessories - Holders and Product Carriers



Holding clamps for Erlenmeyer flasks

- made of stainless steel
- for mounting on the bottom of a Elma cleaning basket
- the Erlenmeyer flask has a firm hold in the holding clamp due to the curved spring



| | 50 ml | 250 ml | 500 ml | 1000 ml |
|--------------------------------|--------------|--------------|----------------|----------------|
| Order number | 1007553 | 1007555 | 1007557 | 1007558 |
| List price | 24,00 € | 30,00 € | 36,00 € | 43,00 € |
| Ext. dimensions W / D / H (mm) | 55 / 55 / 50 | 90 / 90 / 70 | 110 / 110 / 85 | 135 / 135 / 85 |

Dimension definition



All prices ex works, excluding packaging.

The drawings are meant as a visualization aid and are not true to scale.

Technical specifications serve as orientation values and are subject to tolerances due to production and manufacturing. Proper operation can only be guaranteed when using original Elma accessories.

Elmasonic Accessories - Holders and Product Carriers

Holding clamps for Erlenmeyer flasks

Number of Erlenmeyer- / shaking flasks with and without holding clamp in the Elma cleaning baskets

Erlenmeyer- / shaking flasks - without holder

| Passend für Gerätegröße | 30 | 40 | 60 | 80 | 100 | 120 | 150 | 180 | 300 |
|-------------------------|-----------|-----------|-----------|-----------|------------|------------|------------|------------|------------|
| 50 ml (Ø 50 mm) | 6 x | 6 x | 9 x | 16 x | 16 x | 16 x | 40 x | 25 x | 40 x |
| 250 ml (Ø 85 mm) | 2 x | 2 x | 3 x | 5 x | 6 x | 6 x | 15 x | 9 x | 15 x |
| 500 ml (Ø 105 mm) | 1 x | 1 x | 2 x | 4 x | 3 x | 3 x | 8 x | 5 x | 8 x |
| 1000 ml (Ø 131 mm) | -- | -- | -- | -- | 2 x | 2 x | 5 x | 3 x | 5 x |

Erlenmeyer- / shaking flasks - with holder

| Passend für Gerätegröße | 30 | 40 | 60 | 80 | 100 | 120 | 150 | 180 | 300 |
|-------------------------|-----------|-----------|-----------|-----------|------------|------------|------------|------------|------------|
| 50 ml (Ø 55 mm) | 6 x | 5 x | 8 x | 15 x | 14 x | 12 x | 38 x | 20 x | 38 x |
| 250 ml (Ø 90 mm) | 2 x | 2 x | 3 x | 5 x | 5 x | 5 x | 14 x | 8 x | 14 x |
| 500 ml (Ø 110 mm) | -- | -- | 2 x | -- | 3 x | 2 x | 8 x | 5 x | 8 x |
| 1000 ml (Ø 135 mm) | -- | -- | -- | -- | 2 x | 1 x | 5 x | 3 x | 5 x |

Lab bottles (Duran / Schott)

| Passend für Gerätegröße | 30 | 40 | 60 | 80 | 100 | 120 | 150 | 180 | 300 |
|-------------------------|-----------|-----------|-----------|-----------|------------|------------|------------|------------|------------|
| 50 ml (Ø 46 mm) | 6 x | 6 x | 9 x | 16 x | 16 x | 16 x | 40 x | 25 x | 40 x |
| 250 ml (Ø 70 mm) | 2 x | 2 x | 4 x | 7 x | 9 x | 9 x | 18 x | 12 x | 18 x |
| 500 ml (Ø 86 mm) | 2 x | 2 x | 3 x | 5 x | 6 x | 5 x | 14 x | 8 x | 14 x |
| 1000 ml (Ø 101 mm) | 1 x | 1 x | 2 x | 4 x | 4 x | 3 x | 9 x | 5 x | 9 x |

Technical specifications serve as orientation values and are subject to tolerances due to production and manufacturing. Proper operation can only be guaranteed when using original Elma accessories.

Elmasonic Accessories - Holders and Product Carriers

Holding clamps for Erlenmeyer flasks

Example of application



The parts shown are for illustration purposes only and are not included in the scope of delivery.

Technical specifications serve as orientation values and are subject to tolerances due to production and manufacturing. Proper operation can only be guaranteed when using original Elma devices.



Elmasonic Accessories - Holders and Product Carriers

Dental-Tray holder

- made of stainless steel
- Holders for several dental trays 1 / 2 DIN and 3 / 4 DIN

| | 180 |
|--|-----------------|
| Order number | 1104306 |
| List price | 170,00 € |
| Ext. dimensions W / D / H (mm) | 385 / 225 / 165 |
| Int. dimensions W ₁ / D ₁ / H ₁ (mm) | 280 / 205 / 105 |
| Suitable for | |
| • Elmasonic Med | • |

Dimension definition



All prices ex works, excluding packaging.

The drawings are meant as a visualization aid and are not true to scale.

Technical specifications serve as orientation values and are subject to tolerances due to production and manufacturing. Proper operation can only be guaranteed when using original Elma devices.

Elmasonic Accessories - Holders and Product Carriers

Dental-Tray holder

Example of application



The parts shown are for illustration purposes only and are not included in the scope of delivery.

Technical specifications serve as orientation values and are subject to tolerances due to production and manufacturing. Proper operation can only be guaranteed when using original Elma devices.

Elmasonic Zubehör - Inserts

Acid resistant insert tray

- for cleaning with acidic chemicals
- consist of deep-drawn polypropylene and are highly resistant to chemicals
- are hung directly into the ultrasonic unit, which is filled with little water (contact liquid) beforehand.
The ultrasound is then transferred via the contact liquid into the insert tank filled with working liquid (e.g. acid)
- Lid included



| | 30 | 40 | 60 | 100 | 120 | 180 | 300 | 500 |
|--|-----------------|-----------|-----------|-----------------|------------|-----------------|-----------------|-----------------|
| Order number | 2070890000 | | | 2070900000 | | 2070910000 | 2000008000 | 2070780000 |
| List price | 53,00 € | | | 155,00 € | | 205,00 € | 240,00 € | 275,00 € |
| Tank ext. dim. W / D / H (mm) | 300 / 200 / 115 | | | 365 / 305 / 170 | | 390 / 370 / 215 | 570 / 365 / 215 | 530 / 325 / 205 |
| Tank int. dim. W / D / H (mm) | 195 / 95 / 110 | | | 250 / 190 / 165 | | 270 / 245 / 210 | 440 / 240 / 210 | 480 / 275 / 200 |
| Recom. op. capacity V _a (l) | 1,3 | | | 5,8 | | 11,0 | 19,0 | 18,3 |
| Suitable for | 30 | 40 | 60 | 100 | 120 | 180 | 300 | 500 |
| • Elmasonic Easy | • | • | • | • | • | • | • | |
| • Elmasonic Med | • | | • | • | • | • | • | • |
| • Elmasonic Select | • | • | • | • | • | • | • | • |
| • Elmasonic S | • | • | • | • | • | • | • | • |
| • Elmasonic P | • | | • | | • | • | • | |

All prices ex works, excluding packaging.



Technical specifications serve as orientation values and are subject to tolerances due to production and manufacturing. Proper operation can only be guaranteed when using original Elma accessories.

Elmasonic Accessory - Inserts

Acid resistant insert tray



Dimension definition

Tank outer dimensions



Tank internal dimensions

V_a : recommended operating capacity of the cleaning tank
(The actual operating capacity depends on the load.)



The drawings are meant as a visualization aid and are not true to scale.

Technical specifications serve as orientation values and are subject to tolerances due to production and manufacturing. Proper operation can only be guaranteed when using original Elma accessories.

Elmasonic Accessories - Inserts


BD-Steribox

- made of plastic
- with opening in the lid
- with removable sieve for rinsing the contents
- for use in the ultrasonic bath placed directly in the basket or hung over an insert cover



| | |
|--|-----------|
| | |
| Order number | 1105991 |
| List price | 16,50 € |
| Ext. dimensions Dia / H (mm) | 95 / 85 |
| Int. dimensions Dia ₁ / H ₁ (mm) | 75 / 55 |
| Opening dimensions Cover Dia ₃ (mm) | 10 |
| Int. dimensions sieve Dia ₂ / H ₂ (mm) | 70 / 45 |
| Opening dimensions sieve bottom W / D (mm) | 1,8 / 1,8 |
| Capacity max. (ml) | 300 |

Dimension definition



All prices ex works, excluding packaging.

The drawings are merely included as a visualization aid and are not true to scale.

Technical specifications serve as orientation values and are subject to tolerances due to production and manufacturing. Proper operation can only be guaranteed when using original Elma accessories.

Elmasonic Accessories - Inserts

Cleaning cup

- made of plastic
- Lid included
- The cleaning cups can be
 - placed directly in the basket for ultrasonic cleaning devices size 10
 - or hung via an insert cover in the ultrasonic bath.



| | Blue | White |
|--------------------------------|--------------|--------------|
| Order number | 2000001582 | 2000001581 |
| List price | 6,50 € | |
| Ext. dimensions W / D / H (mm) | 95 / 90 / 55 | |
| Diameter Dia (mm) | 75 | |

Dimension definition



All prices ex works, excluding packaging, plus VAT.

The drawings are merely included as a visualization aid and are not true to scale.

Technical specifications serve as orientation values and are subject to tolerances due to production and manufacturing. Proper operation can only be guaranteed when using original Elma accessories.

Elmasonic Accessories - Inserts

Cleaning cup

Example of application



The parts shown are for illustration purposes only and are not included in the scope of delivery.

Technical specifications serve as orientation values and are subject to tolerances due to production and manufacturing. Proper operation can only be guaranteed when using original Elma Units.

Elmasonic Accessories - Inserts

Cleaning glasses


- can be used for simultaneous cleaning of parts with different chemicals, for removing abrasive substances or for cleaning with aggressive chemicals (e.g. acids or alkalis)
- can either be placed directly in a cleaning basket or in combination with an insert cover be hung in the ultrasonic bath
- Plastic lid and rubber ring included




| | 150 ml (without plastic cover) | 600 ml | 600 ml | 1000 ml |
|---|---|-----------------|-----------------|-----------------|
| Order number | 1111834 | 1117632 | 1047272 | 1047376 |
| List price | 15,00 € | 27,00 € | 15,00 € | 38,00 € |
| Glas with cover dimensions ext. Dia / H / D (mm) | 80 / 85 / 95 | 100 / 155 / 110 | 110 / 135 / 125 | 110 / 185 / 125 |
| Glas dimensions internal Dia ₁ / H ₁ (mm) | 55 / 60 | 75 / 145 | 80 / 105 | 85 / 160 |
| Rubber ring external diameter Dia _R (mm) | 70 | | 100 | |

Dimension definition

Cleaning glass with plastic cover



Rubber ring



All prices ex works, excluding packaging.

The drawings are meant as a visualization aid and are not true to scale.

Technical specifications serve as orientation values and are subject to tolerances due to production and manufacturing. Proper operation can only be guaranteed when using original Elma accessories.

Elmasonic Accessories - Inserts

Cleaning glasses

Example of application



The parts shown are for illustration purposes only and are not included in the scope of delivery.

Technical specifications serve as orientation values and are subject to tolerances due to production and manufacturing. Proper operation can only be guaranteed when using original Elma devices.

Elmasonic Accessories - Inserts

Insert cover

- made of stainless steel
- with 2 rubber rings
- to be placed on the ultrasonic cleaning device
- for 2 glasses or 2 cleaning cups



| | 10 | 30 | 40 | 80 | 100 | 120 | 60 | 150 | 180 | 300 | |
|--|----------------|----------------|-----------|-----------|------------|------------|----------------|------------|------------|------------|--|
| Order number | 2070520000 | 2070300000 | | | | | 2070930000 | | | | |
| List price | 30,00 € | 32,00 € | | | | | 38,00 € | | | | |
| External dimensions W / D / H (mm) | 210 / 105 / 10 | 270 / 170 / 15 | | | | | 330 / 180 / 15 | | | | |
| Diameter cover hole Dia (mm) | 85 | 95 | | | | | | | | | |
| External Diameter rubber ring Dia ₁ (mm) | 100 | | | | | | | | | | |
| Suitable for | 10 | 30 | 40 | 80 | 100 | 120 | 60 | 150 | 180 | 300 | |
| • Elmasonic Easy | 1 x | 1 x | 1 x | | 2 x | 2 x | 1 x | | 2 x | 3 x | |
| • Elmasonic Med | | 1 x | | | 2 x | 2 x | 1 x | | 2 x | 3 x | |
| • Elmasonic Select | | 1 x | 1 x | 2 x | 2 x | 2 x | 1 x | 3 x | 2 x | 3 x | |
| • Elmasonic S | 1 x | 1 x | 1 x | 2 x | 2 x | 2 x | 1 x | 3 x | 2 x | 3 x | |
| • Elmasonic P | | 1 x | | | | 2 x | 1 x | | 2 x | 3 x | |


All prices ex works, excluding packaging.

The drawings are meant as a visualization aid and are not true to scale.


Technical specifications serve as orientation values and are subject to tolerances due to production and manufacturing. Proper operation can only be guaranteed when using original Elma accessories.

Dimension definition

Insert cover



Rubber ring




Product Information and Prices

Elma

Elmasonic Accessories - Inserts

Insert cover


- made of stainless steel
- with 2 rubber rings
- for hanging in the ultrasonic bath
- for 2 glasses or 2 cleaning cups




| | 30 | 40 | 60 |
|--|----------------|----|----------------|
| Order number | 1110262 | | 1110261 |
| List price | 50,00 € | | 60,00 € |
| External dimensions W / D / H (mm) | 285 / 110 / 35 | | 350 / 110 / 35 |
| Diameter cover hole Dia (mm) | 65 / 90 | | 95 |
| External Diameter rubber ring Dia ₁ (mm) | 70 / 95 | | 100 |
| Suitable for | 30 | 40 | 60 |
| • Elmasonic Easy | • | • | • |
| • Elmasonic Med | • | | • |
| • Elmasonic Select | • | • | • |
| • Elmasonic S | • | • | • |
| • Elmasonic P | • | | • |

Dimension definition

Insert cover



Rubber ring



All prices ex works, excluding packaging.

The drawings are meant as a visualization aid and are not true to scale.

Technical specifications serve as orientation values and are subject to tolerances due to production and manufacturing. Proper operation can only be guaranteed when using original Elma accessories.

Product Information and Prices



Elmasonic Accessories - Inserts

Insert cover

- made of stainless steel
- with 2 rubber rings
- for hanging in the cleaning basket
- for 2 glasses or 2 cleaning cups



| | 150 | 180 | 300 |
|--|----------------|-----|-----|
| Order number | 1110260 | | |
| List price | 60,00 € | | |
| External dimensions W / D / H (mm) | 260 / 110 / 20 | | |
| Diameter cover hole Dia (mm) | 95 | | |
| External Diameter rubber ring Dia ₁ (mm) | 100 | | |
| Suitable for | 150 | 180 | 300 |
| • Elmasonic Easy | | 2 x | 3 x |
| • Elmasonic Med | | 2 x | 3 x |
| • Elmasonic Select | 3 x | 2 x | 3 x |
| • Elmasonic S | 3 x | 2 x | 3 x |
| • Elmasonic P | | 2 x | 3 x |


All prices ex works, excluding packaging.

The drawings are meant as a visualization aid and are not true to scale.


Technical specifications serve as orientation values and are subject to tolerances due to production and manufacturing. Proper operation can only be guaranteed when using original Elma accessories.

Dimension definition

Insert cover



Rubber ring



Elmasonic Accessories - Inserts

Insert cover


- made of stainless steel
- with 4 rubber rings
- for hanging in the ultrasonic bath
- for 4 glasses or 4 cleaning cups




| | |
|--|----------------|
| | 80 |
| Order number | 1111906 |
| List price | 90,00 € |
| External dimensions W / D / H (mm) | 550 / 130 / 40 |
| Diameter cover hole Dia (mm) | 95 |
| External Diameter rubber ring Dia ₁ (mm) | 100 |
| Suitable for | 80 |
| • Elmasonic Select | • |
| • Elmasonic S | • |

Dimension definition

Insert cover



Rubber ring



All prices ex works, excluding packaging.

The drawings are meant as a visualization aid and are not true to scale.

Technical specifications serve as orientation values and are subject to tolerances due to production and manufacturing. Proper operation can only be guaranteed when using original Elma accessories.


Product Information and Prices



Elmasonic Accessories - Inserts

Insert cover


- made of stainless steel
- with 4 rubber rings
- for hanging in the ultrasonic bath
- for 4 glasses or 4 cleaning cups




| | 100 | 120 | 180 |
|--|----------------|-----|---------------|
| Order number | 1111909 | | 1111921 |
| List price | 60,00 € | | 80,00 € |
| External dimensions W / D / H (mm) | 350 / 230 / 40 | | 380 / 20 / 35 |
| Diameter cover hole Dia (mm) | | 95 | |
| External Diameter rubber ring Dia ₁ (mm) | | 100 | |
| Suitable for | 100 | 120 | 180 |
| • Elmasonic Easy | • | • | • |
| • Elmasonic Med | • | • | • |
| • Elmasonic Select | • | • | • |
| • Elmasonic S | • | • | • |
| • Elmasonic P | • | • | • |

Dimension definition

Insert cover



Rubber ring



All prices ex works, excluding packaging.

The drawings are meant as a visualization aid and are not true to scale.

Technical specifications serve as orientation values and are subject to tolerances due to production and manufacturing. Proper operation can only be guaranteed when using original Elma accessories.

Elmasonic Accessories - Inserts

Insert cover

Example of application



The parts shown are for illustration purposes only and are not included in the scope of delivery.

Technical specifications serve as orientation values and are subject to tolerances due to production and manufacturing. Proper operation can only be guaranteed when using original Elma devices.

Elmasonic Accessories - Inserts

Immersion baskets

- made of plastic and stainless steel
- can - depending on the device size:
 - be hung individually or several times in the ultrasonic bath,
 - be placed individually or several times in a cleaning basket or
 - be immersed individually in a cleaning glass







| Order number | 2070380000 | 1117636 | 2070540000 | 1116797 | 1007545 |
|--|-------------------|------------------|------------------------------|-----------------------------|-------------------|
| List price | 16,00 € | 54,00 € | 51,00 € | 45,00 € | 65,00 € |
| Ext. dim. W / H (mm) | 90 / 180 | 55 / 75 | 75 / 140 | 85 / 105 | 90 / 160 |
| Int. dim. Dia / H ₁ (mm) | 70 / 40 | 45 / 30 | 55 / 60 | 65 / 60 | 70 / 60 |
| Material | plastic | stainless steel | | | |
| Diameter Dia ₁ / Mesh size (mm) | 2,4 | 1x1x0,5 | | | |
| Suitable for cleaning glass | 1000 ml (1047376) | 150 ml (1111834) | 600 ml slim / high (1117632) | 600 ml wide / low (1047272) | 1000 ml (1047376) |

All prices ex works, excluding packaging.

Technical specifications serve as orientation values and are subject to tolerances due to production and manufacturing. Proper operation can only be guaranteed when using original Elma accessories.

Elmasonic Accessories - Inserts

Dimension definition



Example of application



The drawings are meant as a visualization aid and are not true to scale.

The parts shown are for illustration purposes only and are not included in the scope of delivery.

Technical specifications serve as orientation values and are subject to tolerances due to production and manufacturing. Proper operation can only be guaranteed when using original Elma devices.

Product Information and Prices



Elmasonic Accessories - Inserts

Hostal tub

- for cleaning of smallest parts
- protects sensitive parts from damage



| | Small | Big | Big with long temples |
|--|----------------|-----------------|------------------------------|
| Order number | 2070390000 | 207040000 | 2070600000 |
| List price | 19,00 € | 22,00 € | 25,00 € |
| Ext. dimension W / D / H (mm) | 185 / 105 / 95 | 280 / 105 / 100 | 280 / 105 / 125 |
| Int. dimension W ₁ / D ₁ / H ₁ (mm) | 90 / 90 / 55 | | 195 / 90 / 55 |
| Hole diameter Dia (mm) | | 2 | |
| Suitable for | 15 | 20 | 30 |
| • Elmasonic Easy | | • | • |
| • Elmasonic Med | | | • |
| • Elmasonic Select | | | • |
| • Elmasonic S | • | | • |
| • Elmasonic P | | | • |

All prices ex works, excluding packaging.


The drawings are merely included as a visualization aid and are not true to scale.

Technical specifications serve as orientation values and are subject to tolerances due to production and manufacturing. Proper operation can only be guaranteed when using original Elma Units.

Elmasonic Accessories - Inserts

Hostal tub

Dimension definition




The drawings are merely included as a visualization aid and are not true to scale.

Technical specifications serve as orientation values and are subject to tolerances due to production and manufacturing. Proper operation can only be guaranteed when using original Elma Units.

Elmasonic Accessory - Inserts

Silicone mats

- to be placed into the basket for sensitive parts protection
- The silicone mats with order numbers 1082378 and 1117991 must be cut accordingly for the different basket sizes (see cutting guide)
- the silicone mat with order number 1082379 is not cut to size



| | 30 | 60 | 100 | 120 | 150 | 180 | 300 | 300 | 450 | 500 | 600 | 800 | 900 |
|-----------------------------------|---------------|-----------|---------------|------------|------------|------------|------------|------------|------------|------------|---------------|------------|------------|
| Order number | 1082379 | | 1082378 | | | | | | | | 1117991 | | |
| List price | 20,00 € | | 40,00 € | | | | | | | | 108,00 € | | |
| Ext. Dimensions W / D / H (mm) | 193 / 100 / 1 | | 450 / 245 / 1 | | | | | | | | 560 / 460 / 1 | | |
| Mesh size Y x Z (mm) | 5,5 x 1 | | 7,5 x 1 | | | | | | | | 7,5 x 1,5 | | |
| Suitable for | 30 | 60 | 100 | 120 | 150 | 180 | 300 | 300 | 450 | 500 | 600 | 800 | 900 |
| • Elmasonic Easy | • | • | • | • | | • | • | | | | | | |
| • Elmasonic Med | • | • | • | • | | • | • | | • | | | | • |
| • Elmasonic Select | • | • | • | • | • | • | • | • | | • | | | • |
| • Elmasonic S | • | • | • | • | • | • | • | • | | • | | | • |
| • Elmasonic P | • | • | • | • | | • | • | • | | | | | |
| • Elmadry TD | • | | | • | | | • | | • | | • | • | • |
| • Elmasonic xtra ST | | | | | | | | • | | • | • | • | |

All prices ex works, excluding packaging.


Technical specifications serve as orientation values and are subject to tolerances due to production and manufacturing. Proper operation can only be guaranteed when using original Elma accessories.

Elmasonic Accessory - Inserts


Silicone mats

Dimension definition and cutting guide

Order number 1082378



Order number 1117991



The drawings are meant as a visualization aid and are not true to scale.

Technical specifications serve as orientation values and are subject to tolerances due to production and manufacturing. Proper operation can only be guaranteed when using original Elma accessories.

Elmasonic Accessories - Inserts


Small parts holder

- made of stainless steel and Rislan coating to protect delicate parts
- with 10 double hangers, 10 compartments and a subdivided, removable small basket insert (5 compartments à 30x25 mm)



| | 30 | 40 | 60 | 100 | 120 |
|---------------------------------------|-----------|-----------|-----------------|------------|------------|
| Order number | | | 1032059 | | |
| List price | | | 240,00 € | | |
| External dimensions W / D / H (mm) | | | 305 / 115 / 140 | | |
| Height hangers H ₁ (mm) | | | 75 | | |
| Suitable for | 30 | 40 | 60 | 100 | 120 |
| • Elmasonic Easy | 1 x | 1 x | 1 x | 2 x | 2 x |
| • Elmasonic Med | 1 x | | 1 x | 2 x | 2 x |
| • Elmasonic Select | 1 x | 1 x | 1 x | 2 x | 2 x |
| • Elmasonic S | 1 x | 1 x | 1 x | 2 x | 2 x |
| • Elmasonic P | 1 x | | 1 x | 2 x | 2 x |

Dimension definition



All prices ex works, excluding packaging.


The drawings are meant as a visualization aid and are not true to scale.

Technical specifications serve as orientation values and are subject to tolerances due to production and manufacturing. Proper operation can only be guaranteed when using original Elma accessories.

Elmasonic Accessories - Practicals


Cooling coil

- made of stainless steel
- for connection to a closed-circuit recirculating chiller (cryostat)
- for reducing or stabilising the process temperature in the ultrasonic bath
- use of several cooling coils in combination is possible
- optimal for thermostable processes



| | 40 | 70 | 150 | 60 | 80 | 100 | 120 | 130 | 180 | 300 |
|---------------------------------------|----------------|----|-----|----|----|----------------|-----|-----|-----|-----|
| Order number | 1038735 | | | | | 1037515 | | | | |
| List price | 250,00 € | | | | | 250,00 € | | | | |
| External dimensions W / D / H (mm) | 180 / 25 / 130 | | | | | 200 / 25 / 175 | | | | |
| Outer diameter Dia (mm) | | | | 10 | | | | | | |
| Suitable for | 40 | 70 | 150 | 60 | 80 | 100 | 120 | 130 | 180 | 300 |
| • Elmasonic Easy | • | | | | • | | • | | • | • |
| • Elmasonic Med | | | | | • | | • | • | | • |
| • Elmasonic Select | • | | • | • | • | • | • | | • | • |
| • Elmasonic S | • | • | • | • | • | • | • | • | • | • |
| • Elmasonic P | | • | | | • | | | • | • | • |

Dimension definition



All prices ex works, excluding packaging.

The drawings are meant as a visualization aid and are not true to scale.

Technical specifications serve as orientation values and are subject to tolerances due to production and manufacturing. Proper operation can only be guaranteed when using original Elma accessories.

Elmasonic Accessories - Practicals


Drain sieve

- made of stainless steel
- for use in the drain of an Elmasonic ultrasonic device of size 30 to 300
- dirt particles ≥ 0.15 mm are caught in the wire mesh of the drain strainer and thus cannot clog the water drain
- Smallest parts ≥ 0.15 mm are caught with the drain sieve



| Order number | 1045648 |
|------------------------------------|-------------|
| List price | 18,00 € |
| Dimensions ext. Dia / H (mm) | 21 / 35 |
| Diameter sieve Dia_1 (mm) | 10 |
| Mesh size wire mesh (mm) | 0,15 x 0,15 |

Dimension definition



All prices ex works, excluding packaging.

The drawings are meant as a visualization aid and are not true to scale.

Technical specifications serve as orientation values and are subject to tolerances due to production and manufacturing. Proper operation can only be guaranteed when using original Elma accessories.

Elmasonic Accessories - Practicals


Sieve cleaning module SRH 4 / 200

- made of stainless steel and with integrated rotating unit
- up to 4 test sieves [\varnothing up to 200 / 203 mm (8“)] can be cleaned simultaneously in the ultrasonic bath with the SRH 4 / 200. The inclined position of the sieves in the SRH 4 / 200 ensures a better cleaning and rinsing effect
- for gentle cleaning of the test sieves without changing the mesh size or fabric tension. This ensures that the dimensional accuracy and calibration of the sieves is maintained 100 %



| | 300 | 450 | 500 | 900 |
|----------------------|-----|-----|-----------------|-----|
| Order number 230V | | | 1042706 | |
| Order number 115V | | | 1044904 | |
| List price | | | 1.565,00 € | |
| Dimensions W / D / H | | | 675 / 295 / 330 | |
| Weight (kg) | | 6,8 | | |
| Suitable for | 300 | 450 | 500 | 900 |
| • Elmasonic Easy | 1 x | | | |
| • Elmasonic Med | 1 x | | | |
| • Elmasonic Select | 1 x | 1 x | | 2 x |
| • Elmasonic S | 1 x | | 1 x | 2 x |
| • Elmasonic P | 1 x | | | |

Dimension definition



All prices ex works, excluding packaging.

The drawings are meant as a visualization aid and are not true to scale.

Technical specifications serve as orientation values and are subject to tolerances due to production and manufacturing. Proper operation can only be guaranteed when using original Elma devices.

Elmasonic Accessories - Practicals

Sieve cleaning module SRH 4 / 200

Example of application



The parts shown are for illustration purposes only and are not included in the scope of delivery.

Technical specifications serve as orientation values and are subject to tolerances due to production and manufacturing. Proper operation can only be guaranteed when using original Elma devices.

Elmasonic Accessory - Elma Mutebox

Noise protection box

- Noise of the ultrasonic cleaning unit operated in the mute box is reduced by max. 85
- vapors are exhausted via a silenced drain, so that no heat or moisture builds up inside the box.
- build-up of moisture inside the box
- laterally arranged halogen lamps
- Stainless steel jacket
- Product contents:

Elma Mutebox ML: 1x Elma Mutebox ML, 1x Bottom tray with gland plate, 1x Power supply unit,
1x Manual

Elma Mutebox XL: 1x Elma Mutebox XL, 2x Bottom tray with gland plate, 1x splash plate
1x Power supply unit, 1x Manual



Ultrasonic devices not included in the scope of delivery

| Unit size | 10 | 20 | 30 | 40 | 60 | 100 | 120 | 180 | 70 | 80 | 150 | 200 | 300 |
|------------------------|----------|----|----|----|----|-----|-----|-----|----|----|-----------|-----|-----|
| Size | ML | | | | | | | | | | XL | | |
| Order number | 1111901 | | | | | | | | | | 1123006 | | |
| List price | 930,00 € | | | | | | | | | | 1090,00 € | | |
| Suitable for unit size | 10 | 20 | 30 | 40 | 60 | 100 | 120 | 180 | 70 | 80 | 150 | 200 | 300 |
| • Elmasonic Easy | • | • | • | • | • | • | • | • | | | | | • |
| • Elmasonic Med | | | • | | • | • | • | • | | | • | • | |
| • Elmasonic Select | | | • | • | • | • | • | • | | • | • | | • |
| • Elmasonic P | | | • | | • | | • | • | • | | | | • |
| • Elmasonic xtra TT | | | • | | • | | • | | | | | • | |

All prices ex works, excluding packaging.

Elmasonic device in the picture is for illustration only and is not included in the price.

Technical specifications serve as orientation values and are subject to tolerances due to production and manufacturing. Proper operation can only be guaranteed when using original Elma accessories.

Elmasonic Accessory - Elma Mutebox

Noise protection box

Technical Data

| Unit size | 10 | 20 | 30 | 40 | 60 | 100 | 120 | 180 | 70 | 80 | 150 | 200 | 300 |
|--|----|----|----|------------------|----|-----------|-----|-----|------------------|----|-----|-----|-----|
| Ext. Dimensions W / D / H (mm) (Lid and door closed) | | | | 545 / 530 / 525 | | | | | 860 / 605 / 525 | | | | |
| Ext. Dimensions W ₁ / D ₁ / H ₁ (mm) (Lid and door open) | | | | 555 / 975 / 1015 | | | | | 860 / 975 / 1085 | | | | |
| Int. Dimensions W ₂ / D ₂ / H ₂ (mm) | | | | 430 / 380 / 445 | | | | | 710 / 450 / 445 | | | | |
| Weight (kg) | | | | 18,9 | | | | | 26,2 | | | | |
| Mains voltage (V~) | | | | | | 100 - 240 | | | | | | | |

All prices ex works, excluding packaging.

Elmasonic device in the picture is for illustration only and is not included in the price.



Technical specifications serve as orientation values and are subject to tolerances due to production and manufacturing. Proper operation can only be guaranteed when using original Elma accessories.

Elmasonic Accessory - Elma Mutebox



Noise protection box

Dimension definition



External dimensions
(Lid and door closed)



Internal dimensions



External dimensions
(Lid and door opened)



The drawings are meant as a visualization aid and are not true to scale.

Technical specifications serve as orientation values and are subject to tolerances due to production and manufacturing. Proper operation can only be guaranteed when using original Elma accessories.

Elmasonic Accessory - Elma Sink

Rinsing device

- for immersion rinsing or rinsing off objects
- Stainless steel tub
- Elma Sink L and XL are equipped with a drain
- the bath sizes are uniform to Elmasonic Easy, Med, Select and P
- Product contents:
 - Elma Sink S & M: 1x Elma Sink S or M, 1x Operating and installation instructions
 - Elma Sink L & XL: 1x Elma Sink L or XL, 1x Embout de tuyau
1x Collier de serrage, 1x Operating and installation instructions



| Unit size | S | M | L | XL |
|--------------------|---------|----------|----------|----------|
| Order number | 1113839 | 1113840 | 1113841 | 1113843 |
| List price | 90,00 € | 180,00 € | 280,00 € | 340,00 € |
| Tub size | 40 | 60 | 180 | 300 |
| passende Korbgröße | | | | |
| 20 | | 2 x | | |
| 30 | 1 x | | | |
| 40 | 1 x | | | |
| 60 | | 1 x | 2 x | 4 x |
| 80 | | | 1 x | 2 x |
| 100 | | | | 2 x |
| 120 | | | 1 x | 2 x |
| 150 | | | | 1 x |
| 180 | | | 1 x | 1 x |
| 300 | | | | 1 x |

All prices ex works, excluding packaging.

Technical specifications serve as orientation values and are subject to tolerances due to production and manufacturing. Proper operation can only be guaranteed when using original Elma accessories.

Product Information and Prices



Elmasonic Accessory - Elma Sink

Rinsing device

Technical Data - Mechanics

| Unit size | S | M | L | XL |
|--|-----------------|-----------------|-----------------|---|
| Total capacity V (l) | 3,9 | 5,9 | 17,8 | 27,5 |
| Recom. operating capacity V _a (l) | 2,9 | 4,3 | 14,2 | 22,0 |
| Unit ext. dimensions W / D / H (mm) | 305 / 155 / 220 | 365 / 170 / 270 | 395 / 315 / 325 | 565 / 315 / 325 |
| Opening dimensions W ₂ / D ₂ (mm) | 235 / 135 | 300 / 150 | 325 / 295 | 500 / 295 |
| Tank int. dimensions W ₁ / D ₁ / H ₁ (mm) | 215 / 115 / 145 | 280 / 130 / 145 | 300 / 270 / 195 | 475 / 270 / 195 |
| Weight (kg) | 2,0 | 2,8 | 4,7 | 6,1 |
| Connection thread drain (Zoll) | -- | | 3/8 | |
| Drain hose - ID (mm) | -- | | 12 | |
| Basket int. dimensions W ₃ / D ₃ / H ₃ (mm) | 195 / 105 / 50 | 195 / 105 / 95 | 255 / 105 / 95 | 250 / 195 / 90 250 / 195 / 135 280 / 255 / 135 450 / 255 / 55 450 / 255 / 135 |
| Total basket loading (kg) | 1,0 | 3,0 | 5,0 | 6,0 7,0 8,0 7,0 10,0 |



Technical specifications serve as orientation values and are subject to tolerances due to production and manufacturing. Proper operation can only be guaranteed when using original Elma accessories.

Elmasonic Accessory - Elma Sink

Rinsing device



Dimension definition

Unit ext. dimensions




Tank int. dimensions

V_a : recommended operating capacity of the cleaning tank
(The actual operating capacity depends on the load.)



Basket int. dimensions



The drawings are meant as a visualization aid and are not true to scale.

Technical specifications serve as orientation values and are subject to tolerances due to production and manufacturing. Proper operation can only be guaranteed when using original Elma accessories.

Elmasonic Accessory - Elma Dry

Drying units

- from Elma Dry 120 suitable for the use of different basket sizes
- simple and quick operation via a toggle switch with three positions: Off, warm air, cold air
- drying temperature max. 70 °C
- cooling temperature to room temperature
- cover to save energy and to create an effective flow turbulence around the material to be dried
- automatic safety shutdown if temperature is too high
- dust filter (replaceable) prevents dust from the room air from re-contaminating the material to be dried (filter fineness 33ppi)
- stainless steel housing, drip-proof
- very gentle air flow for sensitive parts
- Product contents:
 - 1x Elma Dry 30, 120 or 300, 1x Cover,
 - 1x Power cable with rubber connector (2 m length), 1x Operating manual



| Unit size | 30 | 120 | 300 |
|------------------------------|----------|----------|----------|
| Order number 220 - 240 V~ | 1082175 | 1078542 | 1078575 |
| Order number 115 - 120 V~ | 1082662 | 1082664 | 1082666 |
| List price | 390,00 € | 550,00 € | 990,00 € |

All prices ex works, excluding packaging.

Technical specifications serve as orientation values and are subject to tolerances due to production and manufacturing. Proper operation can only be guaranteed when using original Elma accessories.

Product Information and Prices



Elmasonic Accessory - Elma Dry

Drying units

Technical Data - Mechanics

| Elma Dry | 30 | 120 | 300 |
|--|-----------------|--------------------------|--|
| Unit ext. dimensions W / D / H (mm) | 300 / 186 / 219 | 362 / 285 / 322 | 568 / 347 / 322 |
| Tank int. dimensions W ₁ / D ₁ / H ₁ (mm) | 220 / 120 / 95 | 275 / 215 / 195 | 475 / 270 / 195 |
| Opening dimensions W ₂ / D ₂ (mm) | 235 / 135 | 295 / 235 | 500 / 395 |
| Mesh size Y x Z (mm) | | 2 x 0,63 | |
| Weight incl. cover (kg) | 3,6 | 6,1 | 10,2 |
| Load max. (kg) | 1 | 7 | 12 |
| Suitable basket sizes | | | |
| Elmasonic S, P, Med, Select, Easy | 1 x 30/40 | 2 x 30/40 1 x 100/120 | 3 x 60 1 x 60 + 1 x 100/120 2 x 100/120 1 x 150/300 |
| Elmasonic xtra TT | -- | 2 x 30 1 x 60 | 1 x 30 1 x 60 1 x 120 1 x 200 |
| xtra 30, 50, 70, 150 | -- | 1 x 30 1 x 50 | 3 x 30 1 x 70 1 x 150 |

Technical specifications serve as orientation values and are subject to tolerances due to production and manufacturing. Proper operation can only be guaranteed when using original Elma accessories.

Product Information and Prices



Elmasonic Accessory - Elma Dry

Drying units

Technical Data - Electronics

| Elma Dry | 30 | 120 | 300 |
|------------------------|-----------------------|-----|------|
| Mains voltage (V~) | 220 - 240 / 115 - 120 | | |
| Mains frequency (Hz) | 50 / 60 | | |
| Heating power (W) | 480 | 550 | 1100 |
| Power (W) | 500 | 570 | 1140 |
| Temperature max. (°C) | 70 | | |
| Heating time (min) | ~ 15 | | |
| Air output max. (m³/h) | 160 | 320 | |


Technical specifications serve as orientation values and are subject to tolerances due to production and manufacturing. Proper operation can only be guaranteed when using original Elma accessories.

Elmasonic Accessory - Elma Dry




Drying units

Dimension definition

Unit outer dimensions



Tank internal dimensions



The drawings are meant as a visualization aid and are not true to scale.

Technical specifications serve as orientation values and are subject to tolerances due to production and manufacturing. Proper operation can only be guaranteed when using original Elma accessories.