

# Multispeed

## Polishing unit

### Elma Order Nos

Multispeed	380 030 0000
Grinding wheel	380 021 1000
Additional accessories on request	



Pictured Multispeed

### The large dual spindle polishing machine

Powerful polishing with optimized dust suction.

Further advantages:

- Polishing unit for two brushes
- Two speeds 1400 / 2800 rpm
- Two movable dust protection covers for large pieces
- The polishing area is lit up

Filter pads in thigh drawer for quick and easy replacement. The filter pad can be sent to a refinery to recover the precious metal particles.

### Technical Data

Rotational speed (rpm)	1400/2800
Power (W)	220
Suction power (m <sup>3</sup> /h)	480
Outer dimensions of unit w/d/h (mm / inch)	720 / 350 / 450 / 28 3/10 x 13 4/5 x 17 7/10
Weight (ca. kg / lbs)	38.0 / 83

Technical and optical information are subject to change Elma\_ch\_01/07